

Release 2016.2.7 Enhancement

Laboratory

Laboratory print report format updated

The Laboratory printed report format has been updated to achieve a cleaner, more organized look. Some specifics regarding the changes are listed below.

Interim Report

- The institution name is added to the top of the report, followed by the name of the report.
- The key—Low, High, Critical—prints at the bottom of each page only.
- Laboratory name, address and phone number prints at the bottom of the report.
- Medical director name prints at the bottom of the report.
- The left margin is reduced.
- Current age information is removed.

Cumulative Report

- The phrase Work Copy Only is removed; the Station field is removed
- DOB and Sex data fields are added.
- Laboratory name, address and phone number prints at the bottom of the report.
- Medical director name prints at the bottom of the report.
- The left margin is reduced.
- Test results and reference ranges are closer together.

All Tests by Date - Portrait

- The phrase Work Copy Only is removed; the Station field is removed.
- DOB and Sex data fields are added.
- Laboratory name, address and phone number prints at the bottom of the report.
- Medical director name prints at the bottom of the report.
- The left margin is reduced.
- Test results and reference ranges are closer together.

Lab Status Report

- The institution name is added to the top of the report, followed by the name of the report.
- The phrase **Work Copy Only** is removed; the **Station** field is removed.
- DOB and Sex data fields are added.
- Laboratory name, address and phone number prints at the bottom of the report.
- Medical director name prints at the bottom of the report.

Required action: To enable the new report format, activate a new parameter, MSC LAB ALTERNATE REPORTS, by selecting **Y** for yes at the **Use Alternate Lab Reports** prompt. Test by printing the above lab reports and reviewing for all data pertinent to your normal workflow.

Related artifact: 21125