

Delete an ANATOMIC PATHOLOGY unable to complete

- D Data entry, anat path ...
- E Edit/modify data, anat path ...
- I Inquiries, anat path ...
- L Log-in menu, anat path ...
- P Print, anat path ...
- R SNOMED field references ...
- S Supervisor, anat path ...
- V Verify/release menu, anat path ...
- C Clinician options, anat path ...
- W Workload, anat path ...

Select Anatomic pathology PROD ACCOUNT Option: L Log-in menu, anat path

- LI Log-in, anat path
- DA Delete accession #, anat path
- PB Print log book
- HW Histopathology Worksheet

Select Log-in menu, anat path PROD ACCOUNT Option: LI Log-in, anat path

Select ANATOMIC PATHOLOGY SECTION: MSP SURGICAL PATHOLOGY

Log-In for 2017 ? YES// (YES)

Select Patient Name: TEST,PATIENT 1x-2x-xx 12-xxxx 11xxxxxxxxx Male

<REMOVED PHI>

PATIENT LOCATION: OR// SAMEDAY SURGERY M

Checking surgical record for this patient...

No operations on record in the past 7 days for this patient.

Assign SURGICAL PATHOLOGY (MSP) accession #: 194? YES// N (NO)

Enter Accession #: (1-999999): 1100

SURGEON/PHYSICIAN: USER,PHYSICIAN//

SPECIMEN SUBMITTED BY: USER,PHYSICIAN//

Select SPECIMEN: L4-5,L5-S1//

SPECIMEN: L4-5,L5-S1//

SPECIMEN TOPOGRAPHY: ABDOMEN//

COLLECTION SAMPLE: MISC//

Select SPECIMEN:

DATE/TIME SPECIMEN RECEIVED: SEP 27,2017@14:07//

PATHOLOGIST: USER,PATHOLOGIST//

Select COMMENT:

FROZEN SECTION:

THERE ARE NO LINES!

Edit? NO//

Select LABORATORY TEST: //

This is a required response. Enter '^' to exit

Select LABORATORY TEST: // ^

**No ordered test selected - <ENTRY DELETED>**