

TROUBLE SHOOTING NOTIFICATION ISSUES

OPTION TO GET TO NOTIFICATION MANAGEMENT MENU

Select Systems Manager Menu ED131(GTM) Option: TEST an option not in your menu

Option entry to test: ORB NOT MGR MENU Notification Mgmt Menu

- 1 Enable/Disable Notifications
- 2 Erase Notifications
- 3 Set Urgency for Notifications (GUI)
- 4 Set Deletion Parameters for Notifications
- 5 Set Default Recipient(s) for Notifications
- 6 Set Default Recipient Device(s) for Notifications
- 7 Set Provider Recipients for Notifications
- 8 Flag Orderable Item(s) to Send Notifications
- 9 Archive(delete) after <x> Days
- 10 Forward Notifications ...
- 11 Set Delays for Unverified Orders ...
- 13 Send Flagged Orders Bulletin
- 14 Determine Recipients for a Notification
- 15 Display Patient Alerts and Alert Recipients
- 16 Enable or Disable Notification System
- 17 Display the Notifications a User Can Receive

1) Check to see what users should be getting the notification.

Select Notification Mgmt Menu ED131(GTM) Option: **7 Set Provider Recipients for Notifications**

Set PROVIDER RECIPIENTS Parameters for Notifications

Provider Recipients may be set for the following:

- 1 Division DIV [GENERAL HOSPITAL]
- 2 System SYS [VISTA.GOLD.MEDSPHERE.COM]
- 3 Package PKG [ORDER ENTRY/RESULTS REPORTING]

Enter selection: 3 Package ORDER ENTRY/RESULTS REPORTING

Parameters set for 'Package' may be replaced if ORDER ENTRY/RESULTS REPORTING is installed in this account.

-- Setting Provider Recipients for Package: ORDER ENTRY/RESULTS REPORTING --

Select Notification: ?

Notification	Value
-----	-----
LAB RESULTS	APT
ORDER REQUIRES CHART SIGNATURE	OAPT
FLAG ORDER FOR CLARIFICATION	OAPT
ORDER REQUIRES ELEC SIGNATURE	OAPT

O: ORDERING PROVIDER; A: ATTENDING PHYSICIAN; P: PRIMARY PROVIDER; T: PATIENT CARE TEAMS

2) Check to see how the notification get removed when they are processed

Select Notification Mgmt Menu ED131(GTM) Option: **4 Set Deletion Parameters for Notifications**

Set DELETE MECHANISM Parameters for Notifications

Delete Mechanism may be set for the following:

- 1 Division DIV [GENERAL HOSPITAL]
- 2 System SYS [VISTA.GOLD.MEDSPHERE.COM]
- 3 Package PKG [ORDER ENTRY/RESULTS REPORTING]

Enter selection: 3 Package ORDER ENTRY/RESULTS REPORTING

Parameters set for 'Package' may be replaced if ORDER ENTRY/RESULTS REPORTING is installed in this account.

---- Setting Delete Mechanism for Package: ORDER ENTRY/RESULTS REPORTING ----

Select Notification: ?

Notification	Value
--------------	-------

LAB RESULTS	Individual Recipient
ORDER REQUIRES CHART SIGNATURE	Individual Recipient
FLAG ORDER FOR CLARIFICATION	All Recipients
ORDER REQUIRES ELEC SIGNATURE	All Recipients
ABNORMAL LAB RESULTS (ACTION)	Individual Recipient

INDIVIDUAL RECIPIENT: delete the notification for an individual recipient when a) that individual completes the follow-up action on notifications with associated follow-up action, b) that individual reviews notifications without follow-up actions.

ALL RECIPIENTS: delete the notification for all recipients when a) any recipient completes the follow-up action on notifications with follow-up actions, b) any recipient reviews notifications without follow-up actions.

3) Check which notification a user can Receive

Select Notification Mgmt Menu ED131(GTM) Option: **17 Display the Notifications a User Can Receive**

Enter user's name: VEGA,NICOLE// NV SYSTEM MANAGER

Would you like help understanding the list of notifications? No// (No)

This will take a moment or two, please stand by.....

.....

DEVICE: HOME// TELNET (80x24)

Notification List for VEGA,NICOLE

Page: 1

Notification ON/OFF For This User and Why

ABNL IMAGING RESLT, NEEDS ATTN ON System value is Mandatory
ABNORMAL LAB RESULT (INFO) OFF OERR value is Disabled
ABNORMAL LAB RESULTS (ACTION) OFF OERR value is Disabled
ADMISSION OFF System value is Disabled
ANATOMIC PATHOLOGY RESULTS ON Division GENERAL HOSPITAL value is Mandatory
CCDA TO SEND OFF No Mandatory/Disabled/Enabled values

4) Check Fileman Alert files.

Create a fileman SEARCH to see who received the notification and to check if it was processed and by who.

- Enter or Edit File Entries
- Print File Entries
- Search File Entries
- Modify File Attributes
- Inquire to File Entries
- Utility Functions ...
- Data Dictionary Utilities ...
- Transfer Entries
- Other Options ...
- Electronic Data Dictionary

Select VA FileMan ED131(GTM) Option: **SEARCh File Entries (*MAKE SURE CAPS LOCK IS ON**)**

Output from what File: **ALERT TRACKING//** (28795 entries)

-A- SEARCH FOR ALERT TRACKING FIELD: **PATIENT**
-A- CONDITION: = **EQUALS**
-A- EQUALS PATIENT: TESTPATIENT,ABIGAIL JANE ← **(ENTER PATIENT'S FULL NAME)**

Warning : You have selected a test patient.
WARNING : You may have selected a test patient.

-B- SEARCH FOR ALERT TRACKING FIELD: **DATA FOR PROCESSING**
-B- CONDITION: = **EQUALS**
-B- EQUALS: **10485@** ← **(ENTER THE ORDER NUMBER FOLLOWED BY AN "@")**

-C- SEARCH FOR ALERT TRACKING FIELD:

IF: **AB** PATIENT EQUALS 711 (TESTPATIENT,ABIGAIL JANE)
and DATA FOR PROCESSING EQUALS "10485@"

OR:

STORE RESULTS OF SEARCH IN TEMPLATE:

Sort by: NAME//

Start with NAME: FIRST//

First Print FIELD: **NUMBER**

Then Print FIELD:

Heading (S/C): ALERT TRACKING Search Replace

DEVICE: TELNET (80x24)

ALERT TRACKING Search

OCT 20,2017@09:52 PAGE 1

NUMBER

NAME: **OR,711,12;16;3171016.085908** ← THIS IS THE NAME YOU NEED WHEN DO AN INQUIRY)
30544

1 MATCH FOUND.

DO AN INQUIRY TO THE FILE ALERT TRACKING

Enter or Edit File Entries

Print File Entries

Search File Entries

Modify File Attributes

Inquire to File Entries

Utility Functions ...

Data Dictionary Utilities ...

Transfer Entries

Other Options ...

Electronic Data Dictionary

Select VA FileMan ED131(GTM) Option: **INquire to File Entries**

Output from what File: **ALERT TRACKING//** (28795 entries)

Select ALERT TRACKING NAME: OR,711,12;16;3171016.085908 **OR,711,12;16;3171016.085908**

Another one:

Standard Captioned Output? Yes// (Yes)

Include COMPUTED fields: (N/Y/R/B): NO// **BOTH** Computed Fields and Record Number
(IEN)

NUMBER: 30544

NAME: OR,711,12;16;3171016.085908

DATE CREATED: OCT 16,2017@08:59:08 PKG ID: OR,12

PATIENT: TESTPATIENT,ABIGAIL JANE GENERATED BY: USER,NURSE

GENERATED WHILE QUEUED: YES RETENTION DATE: NOV 30,2017

DISPLAY TEXT: TESTPATIE (T3847): Order requires electronic signature.
ROUTINE TAG: ESORD ROUTINE FOR PROCESSING: ORB3FUP1
DATA FOR PROCESSING: 10485@

****ALL RECIPIENTS RECEIVED THIS NOTIFICATION BECAUSE PROVIDER RECIPIENTS Parameters FOR THIS NOTIFICATION. ONCE THE NOTIFICATION WAS PROCESSED IT DELETED THE NOTIFICATION FROM ALL RECIPIENTS BECAUSE DELETE MECHANISM Parameters WAS SET TO ALL RECIPIENTS FOR THIS NOTIFICATION.*****

RECIPIENT: MANAGER,SYSTEM ← **RECIPIENT WHO RECEIVES THE NOTIFICATION**
RECIPIENT TYPE: INITIAL RECIPIENT ALERT DATE/TIME: OCT 16,2017@08:59:08
RECIPIENT: USER,NURSE DELETED ON: OCT 20,2017@10:09
RECIPIENT TYPE: INITIAL RECIPIENT ALERT DATE/TIME: OCT 16,2017@08:59:08
RECIPIENT: USER,PHYSICIAN DELETED ON: OCT 20,2017@10:09
RECIPIENT TYPE: INITIAL RECIPIENT ALERT DATE/TIME: OCT 16,2017@08:59:08

RECIPIENT: VEGA,NICOLE ← WHO AND WHEN DELETED

ALERT FIRST DISPLAYED: OCT 16,2017@08:59:49

FIRST SELECTED ALERT: OCT 20,2017@10:08:48

PROCESSED ALERT: OCT 20,2017@10:09 DELETED ON: OCT 20,2017@10:09

RECIPIENT TYPE: INITIAL RECIPIENT ALERT DATE/TIME: OCT 16,2017@08:59:08
RECIPIENT: HOOTEN,SARA
ALERT FIRST DISPLAYED: OCT 16,2017@13:48:59
DELETED ON: OCT 20,2017@10:09
RECIPIENT TYPE: INITIAL RECIPIENT ALERT DATE/TIME: OCT 16,2017@08:59:08