

HOW TO SET UP OR CHANGE LAB CHARGE POINTS

- Core Applications ...
- Device Management ...
- Menu Management ...
- Programmer Options ...
- Operations Management ...
- Spool Management ...
- Information Security Officer Menu ...
- Taskman Management ...
- User Management ...
- FM VA FileMan ...
- Application Utilities ...
- Capacity Planning ...
- HL7 Main Menu ...
- Test an option not in your menu

Select Systems Manager Menu PROD ACCOUNT Option: test an option not in your menu
Option entry to test: XPAR EDIT PARAMETER Edit Parameter Values
--- Edit Parameter Values ---

Select PARAMETER DEFINITION NAME: MSC LAB CHARGE POINT

----- Setting MSC LAB CHARGE POINT for System: TESTEMR1.TESTDOMAIN -----
Select LAB CHARGE POINT: ?

LAB CHARGE POINT Value

URINALYSIS ACCESSION
HEMATOLOGY ACCESSION
CHEMISTRY ACCESSION
MICROBIOLOGY ACCESSION

Answer with ACCESSION AREA, or LRDFN, or UID, or HOST UID, or OUTSIDE ACCESSION
Do you want the entire 17-Entry ACCESSION List? N (No)
Select LAB CHARGE POINT: CHEMISTRY

LAB CHARGE POINT: CHEMISTRY// CHEMISTRY CHEMISTRY
Value: ACCESSION// ?

Enter a code from the list.

Select one of the following:

- 0 ACCESSION
- 1 ORDER
- 2 RESULT
- 3 OFF