

Status Time Radiology Report

- Core Applications ...
- Device Management ...
- Menu Management ...
- Programmer Options ...
- Operations Management ...
- Spool Management ...
- Information Security Officer Menu ...
- Taskman Management ...
- User Management ...
- FM VA FileMan ...
- MAIL MailMan Menu ...
- Application Utilities ...
- Capacity Planning ...
- HL7 Main Menu ...
- Test an option not in your menu

Select Systems Manager Menu Option: **TEST an option not in your menu**

Option entry to test: **RAD/NUC MED TOTAL SYSTEM MENU RA OVERALL** Rad/Nuc Med Total System Menu

***** CPT COPYRIGHT NOTICE *****

Any five-digit numeric Physician's Current Procedural Terminology, Fourth Edition (CPT-4) codes, service description, instructions and/or guidelines only are Copyright 1988 (or such other date of publication of CPT-4 as defined in the Berne Implementation Act of 1988 formerly the Copyright Revision Act of 1976) American Medical Association.

No fee schedules, basic unit values, relative value guides or related listings are included in CPT-4. The AMA assumes no responsibility for the consequences attributable to or related to any use or interpretation of information contained in or not contained in this publication. The AMA shall not be deemed to be engaged in the practice of medicine or dispensing medical services.

Printing of any CPT information that will be released external to the VA (excluding areas of billing/fee basis processing, administrative management, clinical management including research, and patient coding/summarizing) must include the following notice: "CPT five-digit codes and/or descriptions only are copyright 1988 AMA (or such other date or publication of the work as defined in the Berne Implementation Act of 1988, formerly the Copyright Revision Act of 1976)."

Press any key to continue

Please select a sign-on Imaging Location: US// GENERAL RADIOLOGY (GE NERAL RADIOLOGY-065)

Default Flash Card Printer: HOME// TELNET (80x40)
Default Jacket Label Printer: HOME// TELNET (80x40)
Default Report Printer: HOME// TELNET (80x40)

Welcome, you are signed on with the following parameters:

Printer Defaults

Version : 5.0 -----
Division : MEMORIAL HOSPITAL ST Flash Card : None
Location : GENERAL RADIOLOGY No cards
Img. Type: GENERAL RADIOLOGY Jacket Label: None
User : VEGA,NICOLE
Report : None

Exam Entry/Edit Menu ...
Films Reporting Menu ...
Management Reports Menu ...
Patient Profile Menu ...
Radiology/Nuclear Med Order Entry Menu ...
Supervisor Menu ...
Switch Locations
Update Patient Record
User Utility Menu ...

Select Rad/Nuc Med Total System Menu Option: **MANAgement Reports Menu**

Daily Management Reports ...
Functional Area Workload Reports ...
Personnel Workload Reports ...
Special Reports ...
Timeliness Reports ...

Select Management Reports Menu Option: **SPEcial Reports**

AMIS Code Dump by Patient
AMIS Report
Camera/Equip/Rm Report
Cost Distribution Report
Detailed Procedure Report
Film Usage Report
Procedure Scaled wRVU/CPT Report
Procedure wRVU/CPT Report

Procedure/CPT Statistics Report
Status Time Report
Wasted Film Report

Select Special Reports Option: **STATUS Time Report**

Select all requesting locations? Y/N: Y

Select **IMAGING TYPE**: GENERAL RADIOLOGY

Select all procedures? Y/N: Y

**** Date Range Selection ****

Beginning DATE : T-7 (JUN 13, 2017)

Ending DATE : T (JUN 20, 2017)

Do you wish to print detailed reports? No//

DEVICE: HOME// TELNET (80x40)

**** Status Tracking Statistics Report **** Page: 1
Division Summary Procedure Detail

Run Date: 06/20/17 For Period: 06/13/17 - 06/20/17