

## Release 2018.1.11 Enhancement

### Radiology

#### MSC Radiology processor now sortable by room/bed

Users can now sort the MSC Radiology processor (MSCRA) by room/bed to facilitate workflow.

```
IM Change Imaging Location          SU Exam Status Update
CA Cancel Exam                     SO Change Sort Order
RG Register Patient for Exam        RF Refresh
LO Change Sign On Imaging Location  DV Detail View
Select Action: Quit// so Change Sort Order
Select sort type: (P/T/E/R/O/B): T// ??

Enter a code from the list.
Select one of the following:

    P      PATIENT
    T      IMAGING TYPE
    E      EXAM STATUS
    R      PROCEDURE
    O      ORDER DATE/TIME
    B      ROOM/BED

Select sort type: (P/T/E/R/O/B): T// █
```

Figure 1: MSC Radiology processor with room/bed sorting option

**Required action:** Test in your normal radiology procedures process.

**Related artifact:** 22324

## Release 2018.1.11 Update

### *Radiology*

#### *MSCRA PROCESSOR more rapidly returns list of patients*

The MSCRA PROCESSOR now returns a list of patients more quickly.

**Required action:** Test the MSCRA PROCESSOR in your normal radiology processing workflow.

**Related artifact:** 22560